

## In situ hybridization for *Gryllus* embryos (10/03/2008)

Fixed embryos are stored in 100% MeOH at -20°C, all steps are done in 1.5 ml tubes  
Prepare PT (2l), MeOH/PT, Hyb solution 1

### Day 1 (5h)

#### Pre-treatment

- 1) **Rehydration**
  - Wash 1x10', 50 % MeOH / PT
  - Wash 2x5', PT
  
- 2) **Permeabilization**
  - Wash 1x 5' in PT-ProK (2µg/ml)
  - Wash 2x 5' in PT
  
- 3) **Post-Fixation**
  - 1x20' n PT-FA
  - Wash 1x 1', 2x 5' in PT

#### Hybridisation

- 4) **Pre-hybridization**
  - 1x 5' PT/50% Hyb solution A
  - 1x 10', 2' 60' (min. 30', can be 2-3h) in 100% hyb solution A at 60°C

Comments:

  - mix embryos with a finger tap a couple of times
  - embryos can be stored in hyb solution at -20C
- 5) **Hybridization**
  - Aliquot embryos to about 50µl per eppendorf
  - Dilute probe in 100ul hyb solution (mostly 1ul of probe stock is sufficient)
  - Boil the probe for 5' to denature; immediately chill on ice
  - Add the probe to the embryos and keep them at 60°C over night (note time )
  - Mix embryos with a finger tap a couple of times.

Prepare for day 2: hyb solution B, PTB

### Day 2

- 6) **Washes**
  - 1x 60 min prewarmed hyb solution B at 60°C
  - 3x 30min Wash 1 (2x SSC) at 60°C
  - 3x 30min Wash 2 (0.2x SSC) at 60°C
  - 3x 10min PT at RT
  - 3x 1h PT at RT
  
- 7) **Blocking step**
  - Wash 1x 10', 1x 60' (or longer) in PTB
  
- 8) **Primary antibody**
  - Incubate PTB/mouse α-DIG-AP (1:1000) O/N at 4°C

### Day 3

9) **Equilibration & AP staining reaction**

- 2x 1 min PT ○ ○
- 6x 30 min PT ○ ○ ○ ○ ○
- prepare NTMT (9ml/probe) and AP buffer
- 1x AP-buffer for 15' ○
- AP-buffer and NBT/BCIP at RT or 4°C in the dark ○

10) **Stopping AP reaction**

- 5x 5min in PT ○ ○ ○ ○ ○

11) **Decolorization (optional)**

- 50% EtOH/PT at RT: variable time ○
- 100% EtOH ○
- 50% EtOH/PT at RT: variable time ○
- 1x 5' in PT ○

Comments: keep examining the embryos using a stereoscope to avoid over-decolorization.  
EtOH reacts with AP reaction solution resulting in precipitation.

12) **Clearization**

- 25% Glycerol/PT until embryos sink ○
- 50% Glycerol ○

Comments: use 50% if embryos should be mounted individually (easier to pipette..)

**Appendix: Solutions**

- PBS(pH7.4)
  - a. 10x PBS 100ml
  - b. DEPC H<sub>2</sub>O 900ml
  
- PT(pH7.4) (PBS + 0.1% Triton X-100)
  - a. PBS 1L
  - b. Triton X-100 1mL

- Keep stirring the prepared solution O/N, because Triton X-100 is viscous.
  
- PT-FA (pH7.4) (PBS + 0.1% Triton X-100 + 4% Formaldehyde)
  - a. PBT 18 ml
  - b. 40% Formaldehyde 2 mL

- Keep stirring the prepared solution O/N, because Triton X-100 is viscous.
  
- PT-ProK (PT+ 2µg/ml Proteinase K)
  - a. ProK stock solution (20mg / mL) 1µL
  - b. PT 10mL

- Prepare for each use and keep on ice. Use stock from Invitrogen
  
- PTB (PT+ 1.5% Roche Blocking reagent)
  - a. Western Blocking Reagent (10%) 1.5 ml
  - b. 10x PBS 1 ml
  - c. 20% Triton X-100 50 µl
  - d. DEPC H<sub>2</sub>O 7.45 ml

- Dissolve and store aliquots at -20°C.
  
- Hybridization Mix A  
 De-ionize Formamide:  
 Add 4g mixed red resin to 40 ml Formamide in 50 ml Falcon tube and shake for 1h, pour formamide in new tube in such a way that red resin stays in old tube

	<b>for 10ml</b>	<b>for 20ml</b>	<b>for 50ml</b>
Roche Blocking solution (final 2%)	2 ml	4 ml	10 ml
Formamide (final 50%)	5 ml	10 ml	25 ml
20xSSC (pH 7.0) (final 5x)	2.5 ml	5 ml	12.5 ml
10% TritonX-100 (final 0.1%)	0.1 ml	0.2 ml	0.5 ml
10% CHAPS (final 0.1%)	0.1 ml	0.2 ml	0.5 ml
0.5M EDTA (final 5mM)	0.1 ml	0.2 ml	0.5 ml
Heparin (50 mg/ml) (final 50 µg/ml)	10 ul	20 ul	50 ul
yeast tRNA (20 mg/ml) (final 1 µg/ml)	0.5 ul	1 ul	2.5 ul
DEPC H <sub>2</sub> O	189.5 ul	379 ul	947.5 ul

Do not vortex yeast tRNA solution.  
 Dissolve and keep them at 60°C until use.

- Hybridization Mix B  
De-ionize Formamide:  
Add 4g mixed red resin to 40 ml Formamide in 50 ml Falcon tube and shake for 1h, pour formamide in new tube in such a way that red resin stays in old tube. SSC might not dissolve completely...

	<b>for 10ml</b>	<b>for 20ml</b>	<b>for 50ml</b>
Formamide (final 50%)	5.0ml	10ml	25ml
20xSSC (pH 7.0) (final 5x)	2.5ml	5.0ml	12.5ml
10% SDS (final 1%)	1.0ml	2.0ml	5.0ml
DEPC H2O	1.5mL	3.0mL	7.5ml

- Wash 1 (2xSSC)  
**for 500mL**  
20xSSC (pH 7.0) (final 2xSSC) 50mL  
10%CHAPS (final 0.1%) 5.0mL  
H2O up to 500mL
- Wash 2 (0.2xSSC)  
**for 500mL**  
20xSSC (pH 7.0) (final 0.2xSSC) 5mL  
10%CHAPS (final 0.1%) 5.0mL  
H2O up to 500mL

- AP-buffer  
**for 10mL**      **for 50mL**  
1M Tris-HCl (pH 9.5) (final: 0.1M) 1.0mL      5.0mL  
5M NaCl (final 0.1M) 0.2mL      1.0mL  
10% TritonX-100 (final: 0.1%) 0.1mL      0.5mL  
1M MgCl<sub>2</sub> (final: 0.05M) 0.5mL      2.5mL  
H2O up to 10mL      50mL

Comments:

- Add 1M MgCl<sub>2</sub> just before use to avoid pH change.

- Staining solution (NBT / BCIP)  
Check stock solution  
**for 10mL**  
AP buffer 10mL  
NBT (50mg/ml) 66µL  
BCIP (50mg/ml) 33µL  
- Protect from light with foil.