

RNA probe synthesis with Roche SP6 or T7 RT enzyme (30-5-08)

Set up a 20 µl reaction with template: 2µl for linearized plasmid, 800ng for PCR product. Keep buffer and reaction at RT during pipetting the reaction.

	total volume	20 µl
H ₂ O (MilliQ is fine)		variable
DNA template (max. 12µl)		variable
Buffer (Roche [Sp6]: 10 X. Roche [T7]: 10 X)		2
NTP labeling mix (either DIG, FITC, BIOTIN)		2
RNase Inhibitor (Roche 03 335 399 001, 40 units/ul)		2
RNA polymerase (T7 or SP6 RT Roche)		2

2 hrs @ 37 °C.

Prepare hyb buffer, RNA gel

DNase treatment

Treat with 1 µl DNase (Megascript) for 15min at 37°C

RNA precipitation

Add 70 µl H₂O (RNase free)

Add 10 µl Na-Acetate 3M pH 5.5 (final 0.3M)

Add 300 µl of ice-cold EtOH (100 %) (3 volumes)

Precipitate for 30' @ -70 °C, spin 30' @ 4°C, 15000 rpm.

Wash pellet with EtOH (70 %, DEPC H₂O), remove residual liquid and air dry (5').

Dissolve pellet in 20 µl DEPC H₂O

Measure concentration

Take out 0.5 µl of the sample and add 3.5µl RNase free H₂O (dilution 1:8)

Spec on nanodrop and run 300ng on gel

(The average total yield of the 20 µl reaction is between 30-40 µg)

Store 10µl of RNA probe at -80 °C

Fill up 9µl with x µl Hyb buffer (depending on concentration, final conc. 100ng/µl) and store at -20 °C.