

dsRNA synthesis for zygotic RNAi**Primers:**

T7	5'-TAATACGACTCACTATAGG-3'
T7-RNAi	5'-gaattgTAATACGACTCACTATAGG-3'
T7-SP6-RNAi	5'-TAATACGACTCACTATAGGattagtgacactataga-3'
T7-T3-RNAi	5'-TAATACGACTCACTATAGGaattaaccctcactaaaggg-3'

Objective

dsRNA is synthesised from a PCR product with T7 promoter sites on each site with the Ambion T7-Megascript Kit with a yield of >6 µg/µl

1. Template DNA with T7 sites on each end**Invitrogen TaqPCRx DNA Polymerase, Recombinant**

total	100.00	Mastermix	1.30
Distilled water to	69.50	µl	90.35
10X PCR Buffer, Minus Mg	10.00	µl 1X	13.00
50 mM MgCl ₂ / H ₂ O	3.00	µl 1.5 mM	3.90
dNTP (2 mM)	10.00	µl 0.2 mM each	13.00
T7-RNAi	3.00	µl 0.3 µM each	3.90
T7-SP6-RNAi	3.00	µl 0.3 µM each	3.90
Template DNA	1.00	10ng/ul	
Taq DNA Polymerase (5 U/µl)	0.50	µl 0.025 units / ul	0.65
		total	128.70 ul
		aliquot	99.00

PCR: 4x 25µl for each synthesis

1x	5'	3'	94°C
35x	30"	94°C	
	30"	58°C	(for T7-RNAi & T7-SP6-RNAi primers)
	1'	72°C	(1min / 1kb length)
1x	5'	72°C	

Combine to 100ul, take out 5ul and run on 1% agarose gel

2. dsRNA synthesis with Ambion Megascript kit

ATP	8	µl
CTP	8	µl
GTP	8	µl
UTP	8	µl
10x buffer	8	µl
PCR template (1ug)	26	µl
H ₂ O (RNase free)	6	µl
enzyme mix	8	µl

3. Use LiCl precipitation (no DNase treatment necessary)

- Add 120µl RNase free H₂O
- Add 120µl LiCl Solution
- Chill at -20C for >30min
- Spin at 13000RPM for 15min
- Wash with 70%EtOH (DEPC H₂O)
- Spin at 13000RPM for 5min
- Dilute in 40µl 1x injection buffer

4. Aligning of dsRNA and 6µg-µl aliquots

- Align dsRNA with heat block from 95-100°C->50°C
- spec and run on gel (1:8 dilution: use 0.5µl dsRNA and add 3.5µl RNase free H₂O)
- fill up to 6 µg/µl
- aliquots at 2µl and store at -20°C